

color. Sulfide itself prevents the reaction if its concentration is very high, in the range of several hundred milligrams per liter. To avoid the possibility of false negative results, use the antimony method to obtain a qualitative result in industrial wastes likely to contain sulfide but showing no color by the methylene blue method. Iodide, which is likely to be present in oil-field wastewaters, may diminish color formation if its concentration exceeds 2 mg/L. Many metals (e.g., Hg, Cd, Cu) form insoluble sulfides and give low recoveries.

Eliminate interferences due to sulfite, thiosulfate, iodide, and many other soluble substances, but not ferrocyanide, by first precipitating ZnS, removing the supernatant, and replacing it with distilled water. Use the same procedure, even when not needed for removal of interferences, to concentrate sulfide. The automated methylene blue method (4500-S²⁻.E) is relatively free from interferences because gas dialysis separates the sulfide from the sample matrix.

1. Apparatus

Glass bottles with stoppers: See 4500-S²⁻.B.1.

2. Reagents

a. Zinc acetate solution: Dissolve 220 g Zn(C₂H₃O₂)₂ · 2H₂O in 870 mL water; this makes 1 L solution.

b. Sodium hydroxide solution (NaOH), 6N.

3. Procedure

a. Put 0.20 mL (4 drops) zinc acetate solution and 0.10 mL (2 drops) 6N NaOH into a 100-mL glass bottle, fill with sample, and add 0.10 mL (2 drops) 6N NaOH solution. Stopper with no air bubbles under stopper and mix by rotating back and forth vigorously about a transverse axis. For the iodometric procedure, use a 500-mL bottle or other convenient size, with proportionally larger volumes of reagents. Vary volume of reagents added according to sample so the resulting precipitate is not excessively bulky and settles readily. Add enough NaOH to raise the pH above 9. Let precipitate settle for 30 min. The treated sample is relatively stable and can be held for several hours. However, if much iron is present, oxidation may be fairly rapid.

b. If the iodometric method is to be used, collect precipitate on a glass fiber filter and continue at once with titration according to the procedure of 4500-S²⁻.F. If the methylene blue method (4500-S²⁻.D) is used, let precipitate settle for 30 min and decant as much supernatant as possible without loss of precipitate. Refill bottle with distilled water, shake to resuspend precipitate, and quickly withdraw a sample. If interfering substances are present in high concentration, settle, decant, and refill a second time. If sulfide concentration is known to be low, add only enough water to bring volume to one-half or one-fifth of original volume. Use this technique for analyzing samples of very low sulfide concentrations. After determining the sulfide concentration colorimetrically, multiply the result by the ratio of final to initial volume. No concentration or pretreatment steps to remove interferences are necessary for 4500-S²⁻.E.

4500-S²⁻ D. Methylene Blue Method

1. Apparatus

a. Matched test tubes, approximately 125 mm long and 15 mm OD.

b. Droppers, delivering 20 drops/mL methylene blue solution. To obtain uniform drops hold dropper in a vertical position and let drops form slowly.

c. If photometric rather than visual color determination will be used, either:

1) *Spectrophotometer,* for use at a wavelength of 664 nm with cells providing light paths of 1 cm and 1 mm, or other path lengths, or

2) *Filter photometer,* with a filter providing maximum transmittance near 660 nm.

2. Reagents

a. Amine-sulfuric acid stock solution: Dissolve 27 g *N,N*-dimethyl-*p*-phenylenediamine oxalate* in an iced mixture of 50 mL conc H₂SO₄ and 20 mL distilled water. Cool and dilute to 100 mL with distilled water. Use fresh oxalate because an old

supply may be oxidized and discolored to a degree that results in interfering colors in the test. Store in a dark glass bottle. When this stock solution is diluted and used in the procedure with a sulfide-free sample, it first will be pink but then should become colorless within 3 min.

b. Amine-sulfuric acid reagent: Dilute 25 mL amine-sulfuric acid stock solution with 975 mL 1 + 1 H₂SO₄. Store in a dark glass bottle.

c. Ferric chloride solution: Dissolve 100 g FeCl₃ · 6H₂O in 40 mL water.

d. Sulfuric acid solution (H₂SO₄), 1 + 1.

e. Diammonium hydrogen phosphate solution: Dissolve 400 g (NH₄)₂HPO₄ in 800 mL distilled water.

f. Methylene blue solution I: Use USP grade dye or one certified by the Biological Stain Commission. The dye content should be reported on the label and should be 84% or more. Dissolve 1.0 g in distilled water and make up to 1 L. This solution will be approximately the correct strength, but because of variation between different lots of dye, standardize against sulfide solutions of known strength and adjust its concentration so 0.05 mL (1 drop) = 1.0 mg sulfide/L.

Standardization—Prepare five known-concentration sulfide standards ranging from 1 to 8 mg/L as described in 4500-S²⁻.A.6, or proceed as follows: Put several grams of clean,

* Eastman Cat. No. 5672 has been found satisfactory for this purpose.

washed crystals of Na₂S · 9H₂O into a small beaker. Add somewhat less than enough water to cover crystals. Stir occasionally for a few minutes, then pour solution into another vessel. This solution reacts slowly with oxygen but the change is insignificant if analysis is performed within a few hours. Prepare solution daily. To 1 L distilled water, add 1 drop of Na₂S solution and mix. Immediately determine sulfide concentration by the methylene blue procedure and by the iodometric procedure. Repeat, using more than 1 drop Na₂S solution or smaller volumes of water, until at least five tests have been made, with a range of sulfide concentrations between 1 and 8 mg/L. Calculate average percent error of the methylene blue result as compared to the iodometric result. If the average error is negative (i.e., methylene blue results are lower than iodometric results), dilute methylene blue solution by the same percentage, so a greater volume will be used in matching colors. If methylene blue results are high, increase solution strength by adding more dye.

g. *Methylene blue solution II*: Dilute 10.00 mL of adjusted methylene blue solution I to 100 mL with reagent water.

3. Procedure

a. *Color development*: Transfer 7.5 mL sample to each of two matched test tubes, using a special wide-tip pipet or filling to marks on test tubes. If sample has been preserved with zinc acetate, shake vigorously before taking subsample. Add to Tube A 0.5 mL amine-sulfuric acid reagent and 0.15 mL (3 drops) FeCl₃ solution. Mix immediately by inverting slowly, only once. (Excessive mixing causes low results by loss of H₂S as a gas before it has had time to react). To Tube B add 0.5 mL 1 + 1 H₂SO₄ and 0.15 mL (3 drops) FeCl₃ solution and mix. The presence of S²⁻ will be indicated by the appearance of blue color in Tube A. Color development usually is complete in about 1 min, but a longer time often is required for fading out of the initial pink color. Wait 3 to 5 min and add 1.6 mL (NH₄)₂HPO₄ solution to each tube. Wait 3 to 15 min and make color comparisons. If zinc acetate was used, wait at least 10 min before making a visual color comparison.

b. *Color determination*:

1) Visual color estimation—Add methylene blue solution I or II, depending on sulfide concentration and desired accuracy,

dropwise, to the second tube, until color matches that developed in first tube. If the concentration exceeds 20 mg/L, repeat test with a portion of sample diluted tenfold.

With methylene blue solution I, adjusted so 0.05 mL (1 drop) = 1.0 mg S²⁻/L when 7.5 mL of sample are used:

$$\text{mg S}^{2-}/\text{L} = \text{no. drops solution I} + 0.1 (\text{no. drops solution II})$$

2) Photometric color measurement—A cell with a light path of 1 cm is suitable for measuring sulfide concentrations from 0.1 to 2.0 mg/L. Use shorter or longer light paths for higher or lower concentrations. This method is suitable for sample concentrations up to 20 mg/L. Zero instrument with a portion of treated sample from Tube B. Prepare calibration curves on basis of colorimetric tests made on Na₂S solutions simultaneously analyzed by the iodometric method, plotting concentration vs. absorbance. A linear relationship between concentration and absorbance can be assumed from 0 to 1.0 mg/L.

Read sulfide concentration from calibration curve.

4. Precision and Bias

In a study by two chemists working in the same laboratory, the standard deviation estimated from 34 sets of duplicate sulfide measurements was 0.04 mg/L for concentrations between 0.2 and 1.5 mg/L. The average recoveries of known additions were 92% for 40 samples containing 0.5 to 1.5 mg/L and 89% for samples containing less than 0.1 mg/L.

5. Quality Control

The quality control practices considered to be an integral part of each method are summarized in Table 4020:I.

6. Bibliography

- POMEROY, R.D. 1936. The determination of sulfides in sewage. *Sewage Works J.* 8:572.
NUSBAUM, I. 1965. Determining sulfides in water and waste water. *Water Sewage Works* 112:113.

4500-S²⁻ E. Gas Dialysis, Automated Methylene Blue Method

1. Apparatus

Automated analytical equipment: An example of the continuous-flow analytical instrument consists of the interchangeable components shown in Figure 4500-S²⁻:2.

The sampler is equipped with a mixer to stir samples before analysis and the gas dialysis membrane, which is maintained at room temperature, separates H₂S from the sample matrix.

2. Reagents

a. *N,N*-dimethyl-*p*-phenylenediamine stock solution: Dissolve 1 g *N,N*-dimethyl-*p*-phenylenediamine dihydrochloride in 500 mL 6*N* HCl. Prepare fresh monthly. Store in an amber bottle.

b. *N,N*-dimethyl-*p*-phenylenediamine working solution: Dilute 190 mL *N,N*-dimethyl-*p*-phenylenediamine stock solution to 1 L. Store in an amber bottle. Prepare weekly.

c. *Ferric chloride stock solution*: Dissolve 13.5 g FeCl₃ · 6H₂O in 500 mL 5*N* HCl. Store in an amber bottle. Prepare fresh monthly.

d. *Working ferric chloride solution*: Dilute 190 mL ferric chloride stock solution to 1 L. Store in an amber bottle. Prepare fresh weekly.

e. *Hydrochloric acid* (HCl), 6*N*.

f. *Sodium hydroxide stock solution* (NaOH), 1*N*.

g. *Sodium hydroxide* (NaOH), 0.01*N*: Dilute 10 mL NaOH stock solution to 1 L.