

50.1.03 - Infant Formulas, Baby Foods, and Enteral Products

**AOAC Official Method 992.06
Vitamin A (Retinol)
in Milk-Based Infant Formula**

**Liquid Chromatographic Method
First Action 1992**

Codex-Adopted-AOAC Method*

(Applicable milk-based infant formulas containing >500 IU vitamin A per reconstituted quart.)

Results of the interlaboratory study supporting the acceptance of the method (milk-based liquid, ready-to-feed):

Mean recovery = 2658 IU vitamin A/L infant formula

$s_r = 129.6$; $s_R = 279.0$; $RSD_r = 4.9\%$; $RSD_R = 10.5\%$

A. Principle

Vitamin A in test portion of infant formula is saponified, partitioned with organic solvent, separated from sample matrix, and quantified by liquid chromatography.

B. Apparatus

(a) *Liquid chromatograph (LC)*.—Capable of pressures up to 20 700 kPa (3000 psi), with injector capable of 100 μ L injections. Operating conditions: eluent flow rate 1.5 \pm 0.2 mL/min; temperature ambient.

(b) *Detector*.—Capable of measuring absorbance at 336 nm, with sensitivity 0.1 AUFS.

(c) *Column*.—4.6 mm id 15 cm stainless steel, packed with 5 μ m silica-based cyano group stationary phase (Sepralyte CN, Analytichem International, Harbor City, CA, USA, is suitable).

(d) *Spectrophotometer*.—Capable of measuring absorbance at 325 nm.

(e) *Shaking water bath*.—Capable of maintaining $70 \pm 2^\circ\text{C}$, variable speed capable of 60 oscillations/min, with test portion area ca 11 14 in. (28 36 cm) (Precision Scientific Model 25 is suitable).

(f) *Glassware*.—(1) 125 mL separatory funnels. (2) 5 mL volumetric flasks. (3) 100 mL low-actinic volumetric flasks.

C. Reagents

(a) *Mobile phase solution*.—Hexane-isopropyl alcohol (100 + 0.25, v/v), HPLC grade solvents. Degas 2-5 min under vacuum.

(b) *Wash solution*.—H₂O-absolute ethanol (3 + 2, v/v).

(c) *Extraction solution*.—Hexane-methylene chloride (3 + 1, v/v), HPLC grade solvents.

(d) *Saponification solution*.—10.5M potassium hydroxide (KOH). Dissolve 673 g KOH in 1 L H₂O.

(e) *Antioxidant solution*.—1% pyrogallol. Dissolve 5.0 g pyrogallol (1,3,5-trihydroxybenzene, 98%, Aldrich is suitable source) in 500 mL absolute ethanol.

(f) *Standard solutions*.—(1) *Stock standard solution*.—10 mg/mL retinyl palmitate in hexane. Quantitatively transfer equivalent of 0.1 g retinyl palmitate with hexane (HPLC grade) into 100 mL low-actinic volumetric flask, dilute to volume with hexane, and shake well to dissolve. Make fresh every 2 weeks. Store at -20°C in explosion-proof freezer when not in use.

(2) *Intermediate standard solution*.—Pipet 2 mL of stock standard solution, (1), into 250 mL volumetric flask and dilute to volume with hexane.

(3) *Working standard solution*.—Approximately 1.6 $\mu\text{g/mL}$ retinyl palmitate. Pipet 2 mL intermediate standard solution, (2), into 100 mL low-actinic volumetric flask. Evaporate to dryness under nitrogen. Dissolve residue in antioxidant solution, C(e), and dilute to volume. Prepare fresh daily.

D. Extraction of Standards and Test Portions

Pipet 10.0 mL working standard solution, **C(f)(3)**, or test portion containing ca 20 IU vitamin A activity (10 mL for ready-to-feed formulas) into 150 mL centrifuge tube. Bring test portion to 10 mL with H₂O, if necessary. To standard tubes, add 10 mL H₂O, 20 mL antioxidant solution, **C(e)**, and 5 mL saponification solution, **C(d)**. To test portion tubes, add 30 mL antioxidant solution and 5 mL saponification solution. Cap tubes and swirl briefly to mix. Place tubes in 70°C shaking H₂O bath (ca 60 oscillations/min) for 25 min. Remove tubes and place in ice 5 min, or until contents cool to room temperature.

Quantitatively transfer contents to separate 125 mL separatory funnels. Wash remaining test portion or standard from tube into funnel with 30 mL H₂O. Pipet 30.0 mL extraction solvent, **C(c)**, into funnel and shake ca 2 min. When layers separate, discard aqueous (lower) layer. Add 30 mL wash solution, **C(b)**, to funnel and shake very gently 30 s, venting frequently. Let phases separate and discard aqueous layer. Repeat wash step 3. Pipet 20.0 mL portion from funnel to 50 mL tube and evaporate to dryness under nitrogen. Transfer residues quantitatively to separate 5 mL volumetric flasks and dilute to volume with mobile phase solution, **C(a)**.

E. Determination of Standard Concentration

Pipet 2 mL of intermediate standard solution, **C(f)(2)**, into 50 mL volumetric flask and dilute to volume with hexane. Transfer portion of this solution into 1 cm cell path length cuvet and measure absorbance at 325 nm. Calculate concentration of working standard solution, C_{std} , as follows:

$$C_{\text{std}} =$$

where A_{325} = absorbance of working standard solution at 325 nm; $\epsilon = 975$, extinction coefficient of retinyl palmitate in hexane at 325 nm; and $b = 1$ cm, cell path length.

F. System Suitability Test

Inject 100 μ L of saponified working standard solution into LC. Typical retention times for 13-*cis*-retinol and *trans*-retinol are 7.5 and 9.0 min, respectively. Calculate R factor between 13-*cis*-retinol and *trans*-retinol as follows:

$$R =$$

where t_1 and t_2 = retention time measured from injection time to elution time of peak maximum of 13-*cis*-retinol and *trans*-retinol, respectively, and W_1 and W_2 = width of peak measured by extrapolating relatively straight sides to baseline of 13-*cis*-retinol and *trans*-retinol, respectively.

If R factor is <1.3 , increase amount of isopropyl alcohol added per liter [mobile phase solution, **C(a)**] by ca 0.05%. Inject saponified working standard solution 5. Calculate reproducibility of replicate injections in terms of standard deviations (per USP), which should be $\leq 2\%$. Typical relative standard deviation values for peak height are $\pm 1.3\%$.

G. Liquid Chromatography

Inject 100 μ L standard or test solution into LC.

H. Calculations

Since 13-*cis* vitamin A palmitate is not readily available, standard curve for all-*trans* vitamin A palmitate is used to determine biological potencies for both, correcting for 13-*cis* vitamin A palmitate, at 0.75 potency relative to all-*trans* vitamin A palmitate. This is based on assumption that relative molar absorptivities of both isomers are virtually equal at 336 nm.

Measure peak areas of *cis*- and *trans*-isomers of retinol in both test solution and standard chromatograms. Multiply peak area under 13-*cis* vitamin A palmitate curve by 0.75. Sum the 2 areas to represent total peak area. Calculate IU per reconstituted quart of vitamin A activity (V) as follows:

$$V =$$

where A_{sample} = total peak area of test solution; A_{standard} = total peak area of standard; C_{standard} = concentration of working standard solution, $\mu\text{g/mL}$; $1/0.55$ = IU/ μg retinyl palmitate; and 946.33 = mL/quart.

Reference:

J. AOAC Int. **76**, 399(1993).

CAS-68-26-8 (vitamin A)

Revised: March 1997

* Adopted as a Codex Reference Method (Type II) for liquid chromatography of vitamin A (retinol) in infant formula and follow-up formula.

© 2002 AOAC INTERNATIONAL