
50.1.04 - Infant Formulas, Baby Foods, and Enteral Products

**AOAC Official Method 992.03
Vitamin E Activity (All-*rac*- -Tocopherol)
in Milk-Based Infant Formula**

**Liquid Chromatographic Method
First Action 1992
Final Action 1996**

Codex-Adopted-AOAC Method*

(Applicable to determination of vitamin E activity in milk-based infant formula.)

Results of the interlaboratory study supporting the acceptance of the method (milk-based liquid, ready-to-feed):

Mean recovery = 24.14 IU vitamin E/L infant formula

$s_r = 2.04$; $s_R = 2.82$; $RSD_r = 8.46\%$; $RSD_R = 11.69\%$

A. Principle

Vitamin E activity in test portions of infant formula is determined by saponification of all-*rac*- -tocopherol, partitioning with organic solvent, separation from product matrix, and quantification by liquid chromatography.

B. Apparatus

(a) *Liquid chromatograph (LC)*.—Capable of pressures up to 3000 psi with injector capable of 100 μ L injections. Operating conditions: eluent flow rate 2.0 ± 0.2 mL/min; temperature ambient.

(b) *Detector*.—Capable of measuring absorbance at 280 nm, with sensitivity 0.02 AUFS.

(c) *Precolumn*.—2 mm id \times 2 cm stainless steel, packed with 40 μ m pellicular reversed-phase C₁₈ (Alltech 28551 is suitable).

(d) *Column*.—4.6 mm id 25 cm stainless steel, packed with 5 μ m silica (Hypersil Silica is suitable).

(e) *Shaking water bath*.—Capable of maintaining $70 \pm 2^\circ\text{C}$, with variable speed capable of 60 oscillations/min, ca 11 14 in. (28 36 cm) test sample area (Precision Scientific Model 25 is suitable).

(f) *Glassware*.—(1) 125 mL separatory funnels. (2) 5 mL volumetric flasks. (3) 100 mL low-actinic volumetric flasks.

C. Reagents

(a) *Mobile phase solution*.—Hexane-isopropyl alcohol (99.92 + 0.08, v/v), HPLC grade solvents. Degas 2-5 min under vacuum.

(b) *Wash solution*.—H₂O-absolute ethanol (3 + 2, v/v).

(c) *Extraction solution*.—Hexane-methylene chloride (3 + 1, v/v), HPLC grade solvents.

(d) *Saponification solution*.—10.5M potassium hydroxide (KOH). Dissolve 673 g KOH in 1 L H₂O.

(e) *Antioxidant solution*.—1% pyrogallol. Dissolve 5.0 g pyrogallol (1,3,5-trihydroxybenzene, 98%; Aldrich is suitable source) in 500 mL absolute ethanol.

(f) *Standard solutions*.—(1) *Stock standard solution*.—0.5 mg/mL all-*rac*- α -tocopheryl acetate in hexane. Accurately weigh ca 50 mg all-*rac*- α -tocopheryl acetate (USP Reference Standard) into 100 mL low-actinic volumetric flask and dilute to volume with hexane (HPLC grade). Shake well to dissolve. Make fresh every 3 weeks. Store at -20°C in explosion-proof freezer when not in use.

(2) *Working standard solution*.—10 μ g/mL all-*rac*- α -tocopheryl acetate. Pipet 2 mL stock standard solution, (1), into 100 mL low-actinic volumetric flask. Evaporate to dryness under nitrogen. Dissolve residue in antioxidant solution, C(e), and dilute to volume. Prepare fresh daily.

(g) *Suitability test solution*.—Approximately 15 μ g/mL all-*rac*- α -tocopherol (USP Reference Standard) and all-*rac*- α -tocopheryl acetate (USP Reference Standard) in hexane (HPLC grade).

D. Extraction of Standard and Test Portions

Pipet 10.0 mL working standard solution, **C(f)(2)**, or test portion containing ca 0.095 IU vitamin E activity (10 mL for ready-to-feed formulas) into 150 mL centrifuge tube. Bring test solution volume to 10 mL with H₂O, if necessary. To standard tubes, add 10 mL H₂O, 20 mL antioxidant solution, **C(e)**, and 5 mL saponification solution, **C(d)**. To sample tubes, add 30 mL antioxidant solution and 5 mL saponification solution. Cap tubes and swirl briefly to mix. Place tubes in 70°C shaking H₂O bath (ca 60 oscillations/min) for 25 min. Remove tubes and place in ice 5 min, or until contents cool to room temperature.

Quantitatively transfer contents to separate 125 mL separatory funnels. Wash remaining test portion solution or standard from tube into funnel with 30 mL H₂O. Pipet 30.0 mL extraction solvent, **C(c)**, into funnel and shake ca 2 min. When layers separate, discard aqueous (lower) layer. Add 30 mL wash solution, **C(b)**, to funnel and shake very gently 30 s, venting frequently. Let phases separate and discard aqueous layer. Repeat wash step 3. Pipet 20.0 mL portion from funnel into 50 mL tube and evaporate to dryness under nitrogen. Transfer residues quantitatively to separate 5 mL volumetric flasks and dilute to volume with mobile phase solution, **C(a)**. Inject 100 µL standard or sample into LC.

E. System Suitability Test

Inject 100 µL suitability test solution, **C(g)**, into LC. Typical peak retention times for tocopherol and tocopheryl acetate are 6.0 and 5.0 min, respectively. Calculate resolution (*R*) factor between tocopherol and tocopheryl acetate as follows:

$$R =$$

where t_1 and t_2 = retention time measured from injection time to elution time of peak maximum of tocopherol and tocopheryl acetate, respectively, and W_1 and W_2 = width of peak measured by extrapolating relatively straight sides to baseline of alcohol and acetate, respectively.

If *R* factor is >1.0, proceed with test portion analysis; if *R* factor is <1.0, decrease amount of isopropyl alcohol added per liter [mobile phase solution, **C(a)**] by ca 0.01%. Inject working standard solution, **C(f)(2)**, 5. Calculate reproducibility of replicate injections in terms of standard deviations (per USP), which should be ≤ 2%. Typical relative standard deviation values for peak height are ±1.5%.

F. Liquid Chromatography

Inject 100 µL standard or test solution into LC.

G. Calculations

Measure peak heights or peak areas of all-*rac*- α -tocopherol in test solution and standard chromatograms. Calculate IU per reconstituted quart of vitamin E activity (A) as follows:

$$A = H_{\text{sam}}/H_{\text{std}} \quad C_{\text{std}} \quad 0.001 \text{ IU}/\mu\text{g} \quad 946.33 \text{ mL/quart}$$

where H_{sam} = peak height of test portion solution; H_{std} = peak height of standard; C_{std} = concentration of standard, $\mu\text{g/mL}$.

Reference:

J. AOAC Int. **76**, 398(1993).

CAS-59-02-9 (α -tocopherol)

CAS-7695-91-2 (α -tocopheryl acetate)

Revised: March 1997

* Adopted as a Codex Reference Method (Type II) for liquid chromatography of vitamin E (milk based infant formula) in special foods.