

---

**44.5.14 - Sugars and Sugar Products / Maple, Sap, Maple Syrup, Maple Syrup Products**

## **AOAC Official Method 964.21 Formaldehyde in Maple Syrup**

### **Spectrophotometric Method First Action 1964 Final Action 1965**

(Not suitable for beet or cane sugars.)

#### ***A. Apparatus***

(a) *Distillation apparatus*.—30 mL micro Kjeldahl flask, with standard taper 19/38 joint and 10 cm H<sub>2</sub>O-cooled West condenser with standard taper 19/38 inner joint bent at 90° angle. Insulate condenser from standard taper joint to condenser H<sub>2</sub>O jacket. Install adjustable flask heater equipped with conical porcelain coil element and hand-controllable rheostat placed on underside (Cenco, Inc., 2600 Kostner Ave, Chicago, IL 60623, USA, No. 16530-1, or equivalent).

(b) *Spectrophotometer*.—Suitable for measuring *A* at 415 nm; with matched 1 cm cells or matched test tubes.

#### ***B. Reagents***

(a) *Nash's Reagent "B"*.—Dissolve 150 g NH<sub>4</sub>CH<sub>3</sub>COO, 3 mL CH<sub>3</sub>COOH, and 2 mL acetylacetone in 200-300 mL H<sub>2</sub>O in 1 L flask and dilute to volume.

(b) *Formaldehyde*.—Approximately 37% by weight. Assay by [898.01B](#) (see 7.8.04).

#### ***C. Preparation of Standard Solutions***

(a) *Stock solution*.—1000 µg/g. Weigh 5.35 g 37.4% CH<sub>2</sub>O solution (for other concentrates, g CH<sub>2</sub>O solution required = 200/% CH<sub>2</sub>O) into 2 L volumetric flask containing H<sub>2</sub>O and dilute to volume with H<sub>2</sub>O.

(b) *Solution B*.—50  $\mu\text{g/g}$ . Pipet 10 mL stock solution into 200 mL volumetric flask and dilute to volume with  $\text{H}_2\text{O}$ .

(c) *Solution C*.—100  $\mu\text{g/g}$ . Pipet 10 mL stock solution into 100 mL volumetric flask and dilute to volume with  $\text{H}_2\text{O}$ .

(d) *Solution D*.—200  $\mu\text{g/g}$ . Pipet 10 mL stock solution into 50 mL volumetric flask and dilute to volume with  $\text{H}_2\text{O}$ .

(e) *Formaldehyde standard solutions*.—Prepare 1, 2, and 4  $\mu\text{g/g}$  standard solutions by pipetting 10 mL *Solutions B, C, and D*, respectively, into 500 g syrup and stirring mechanically 15 min.

#### ***D. Determination***

Weigh  $20 \pm 0.20$  g syrup into tared 30 mL micro Kjeldahl flask and insulate flask neck with asbestos cord as in **A(a)** from base of standard taper joint down to body of flask. (This insulation is temporary; repeat after each washing and weighing.) Add 2 drops antifoam agent, and connect West condenser. Mount apparatus, adjusting slope of condenser at  $45^\circ$  angle, with flask bottom centered in conical cavity and at such height that syrup level in flask is even with top plate of heater, without touching hot element at any point; *see* Figure **964.21**. Heat flask with heater, previously adjusted with control knob in such position that exactly 3 mL  $\text{H}_2\text{O}$  is distilled from syrup in 12-14 min. Collect 3 mL distillate in 5 mL graduate with funnel top. Using transfer pipets, place 1 mL distillate in 13 mm id test tube, and add 1 mL  $\text{H}_2\text{O}$  and 2 mL Nash's reagent. Heat 30 min in  $\text{H}_2\text{O}$  bath at  $37 \pm 1^\circ\text{C}$  to develop color. Transfer colored solution to 1 cm cell and measure  $A$  at 415 nm against  $\text{H}_2\text{O}$ .

#### ***E. Blanks***

To determine  $A$  due to reagents, substitute 1 mL  $\text{H}_2\text{O}$  (from same source as used in determination) for 1 mL test distillate. Subtract  $A$  of blank from that for sample to obtain  $A$  due to  $\text{CH}_2\text{O}$ . Or, as simpler operation, measure  $A$  of distillate with instrument adjusted to 0  $A$  for blank.

Obtain concentration of  $\text{CH}_2\text{O}$  in syrup from  $A$ , using standard curve.

#### ***F. Preparation of Standard Curve***

Construct standard curve by plotting  $A$  obtained for syrups containing 1, 2, and 4 ppm  $\text{CH}_2\text{O}$  against concentration of added  $\text{CH}_2\text{O}$  in ppm.

Straight line relationship is obtained for standard curve. Project this line to *Y* axis (*A*); *Y* intercept indicates blank for syrup. Since syrup used to construct curve from *A* values may be atypical, draw and use parallel curve with 0 intercept. Correct ppm values obtained from this curve for average syrup blank (ca 0.9 ppm).

**References:**

*JAOAC* **47**, 548(1964); **56**, 132(1973).

CAS-50-00-0 (formaldehyde)